

TRANSTEK

# Cellular BB1597-G

Compact 4G BPM with triple mode



# Cellular Connectivity

 Smart retransmission


 Intelligent Guidance

Pls go to a good signal



# Memory Mode

Memory Mode can store 1×60 measurements, making blood pressure trends visible for easier data management and physician assessment.



User			
2022 06/14	SYS 112	DIA 66	PUL 63
2022 06/15	SYS 111	DIA 62	PUL 65
2022 06/16	SYS 127	DIA 76	PUL 72
2022 06/17	SYS 116	DIA 68	PUL 63



# Medical-grade Accuracy




Clinically Validated Accuracy

Triple Mode Measurement

Intelligent Filtering

# Easy to Read

BP Classification Indicator



LCD VA:




73mm

49mm

Blood pressure classification  
reference: AHA 2025

# Device Management

- Encrypted Data Transmission
- SSL Two-way Digital Certificate+API Key
- Update the Configuration Remotely
- BPM-Fully Devices & BPM-Unmanaged Devices



# Product Parameters

TRANSTEK

Model Number	BB1597-G
Type	Upper Arm
Dimensions (mm)	118×126×72mm
Cuff Size (cm)	22-45cm
Display VA (mm)	73×49mm
Display Type	LCD/RD3 negative
Backlight Type	White Backlight
Wireless	CatM:B2/B4/B12/B13/B25
BP range	SYS 60mmHg ~ 230mmHg ; DIA 40mmHg ~ 130mmHg
Heart rate range	40 ~ 199
Memory Capacity	1*60
Power supply	4 × AA or 6V@1A
Battery life times	≥400
Cuff life times	≥10000
Device life times	≥15000 / 5 years
FDA	K241351

# TRANSTEK

## Your Reliable BP Monitor

[inquiry@transtekcorp.com](mailto:inquiry@transtekcorp.com) | [www.transtekcorp.com](http://www.transtekcorp.com)