

TeleRPM™

TeleRPM BGM(BLE) 2022

Monitor blood glucose easily



8 Electrodes
Auto Coding



One-click
Operation



Rechargeable
Battery



Medical-grade
Accuracy



Product Name	TeleRPM BGM(BLE) 2022
Dimensions	83mm x 52mm x 18.7mm
Display Size	32mm x 32mm
Weight	Approximately 53g
Sample	Fresh capillary whole blood
Test Measured	Glucose in fingertip capillary whole blood
Enzyme	GOD/Glucose Oxidase
Sample Volume	0.8 µL
Test Time	5 seconds
Memory	Up to 500 records with date and time
Measurement range	20-600 mg/dL
Battery Type	Rechargeable, non-serviceable, 250mAh, 3.7 Volt DC nominal, lithium polymer battery (5V input charge voltage)
Automatic Shutoff	2 minutes after last action
Operating temperature	41-113°F
Operating relative humidity	10-90%
Data Port	Micro USB



Contact us

inquiry@transtekcorp.com
www.transtekcorp.com

LinkedIn



Follow us on : Transtek Medical