

# TeleRPM™

## TeleRPM BPM(BLE) 2022

Take charge of your blood pressure easily



9"-17" Arm  
Circumference



One-click  
Operation



Backlit  
LCD Display



Bluetooth 5.0



Product Name	TeleRPM BPM(BLE) 2022
Dimensions	100 x 132 x 48mm
Cuff Size	22-42cm
LCD Display	61 x 74mm
Gross Weight	675g
Measurement Method	Oscillographic measurement
Rated Cuff Pressure	0~299 mmHg(0.0~39.9kPa)
Measurement Range	SYS: 60 ~ 230 mmHg(8.0~30.7kPa) DIA: 40 ~ 130 mmHg(5.3~17.3kPa) Pulse value: 40~199 beats/min
Accuracy	Pressure: $\pm 3$ mmHg, Pulse: $\pm 5\%$
Normal Working Condition	Temperature: +20 °C to +60 °C Relative humidity: 15% ~ 90% Pressure: 700hPa ~ 1060hPa
Storage & Transportation Condition	Temperature: -20 °C ~ 60 °C, Relative humidity: $\leq 93\%$ RH
Power Supply	4 x AAA alkaline batteries, Micro-USB cable
Ingress Protection	IP21



LinkedIn

## Contact us

[inquiry@transtekcorp.com](mailto:inquiry@transtekcorp.com)  
[www.transtekcorp.com](http://www.transtekcorp.com)



Follow us on : Transtek Medical